




Osteogenesis and angiogenesis of a bulk metallic glass for biomedical implants

Sun et al. Bioactive Materials 2022


Allan DIZET & Matteo MARENGO
3A BIOMED - PHELMA
15/12/2022

Implants in orthopedic surgery



← Bone fracture fixation

Dental implant system →



What is important ?


1) Biocompatibility

2) Mechanical behavior

3) Long-term safety / stability

Biocompatibility

What is the reaction when implanting a foreign-body ?



Vroman effect

Protein with the highest mobility will adsorb first.


Bulk Metallic Glasses (BMGs)




A new type of material : BMGs

- Co-Cr alloys / 316L stainless steel
- Toxic effect

AMORPHOUS ATOMIC
STRUCTURE



CRYSTALLINE ATOMIC
STRUCTURE



- Amorphous solid material
- Glass-like atomic structure
- Irregular, very strong and elastic

The three protagonists

Ti-based BMGs



+ Improved wear resistance

- Young's modulus 80-120 GPa > human bone 10-40 GPa



$\text{Zr}_{16}\text{Ti}_2\text{Cu}_{25}\text{Al}_{12}$ (ZT1)

(poly-ether-ether-ketone) PEEK

Cp-Ti (pure titanium)



What is angiogenesis and osteogenesis ?

Osteogenesis and angiogenesis

Osteogenesis : formation of bone

Angiogenesis : formation of new blood vessels from pre-existing vessels




Therefore, what is the problematic of this study ?

Problematic


*“ Is **ZT1** a good candidate for orthopedic applications regarding its biocompatibility and biomechanical properties in in-vitro and in-vivo measurements ? ”*

Challenges of the study

- 1 Characterization of the materials**
- 2 Osteogenesis properties**
- 3 Angiogenesis properties**




Characterization of the materials




In-vitro measurements



In-vivo measurements



Discussion




Characterization of the materials



In-vitro measurements




In-vivo measurements




Discussion

Materials preparation




Rods preparation
→ *in-vivo* evaluations



Plates preparation
→ hydrophilic characterization
→ bio-corrosion tests
→ *in-vitro* evaluations





Structure study

- Glassy structure characterized by X-ray diffraction (XRD)
→ typical broad scattering signal




Amorphous structure


Hydrophilic characterization



ZT1 displays the most hydrophilic surface characteristic




Characterization of the materials



In-vitro measurements



In-vivo measurements



Discussion

What is Cytocompatibility ?

- **cytocompatibility** : property of not being harmful to the cell

How to assess Cytocompatibility ?

- Cell viability
- Cell adhesion
- Cell morphology
- Cell proliferation



Study of cytocompatibility with MC3T3-E1 pre-osteoblasts


MC3T3-E1
pre-osteoblasts



Seeding
Incubation



MC3T3-E1
pre-osteoblasts





What is MC3T3-E1 pre-osteoblasts ?

- **MC3T3** = Osteoblast precursor cell line (from mouse skull)

- **MC3T3-E1 sub-line** = spontaneously transformed (immortalised) cell line.

- Caution when extrapolating these results to normal cells

Cell viability





- live/dead staining assay
- living cells in green
- dead cells in red



ZT1 has the better cytocompatibility

Cell proliferation





- cell proliferations after 2 days / 3 days in CCK-8 assay.
- the darker / the more it proliferates

→ **ZT1 stimulates proliferation**



Angiogenesis measurements with HUVECs







What is HUVECs ?

- Human Umbilical Vein Endothelial cells

- Model system for the study of the function and pathology of endothelial cells (e.g. angiogenesis)

- Easily to proliferate
- Low cost

Angiogenesis assessment




- wound narrowing after 24h for ZT1
- higher transwell migration for ZT1
- tube formation for ZT1


ZT1 enhances angiogenic ability

ZT1 has better **cytocompatibility** and **angiogenesis** than CP-Ti and PEEK

 *In-vivo* measurements




Characterization of the materials



In-vitro measurements




In-vivo measurements





Discussion


- 1 Methods & Experiments**
- 2 Pain measurements**
- 3 Osteogenesis measurements**
- 4 Angiogenesis measurements**




Female Sprague-Dawley rats











Where is the implant ?



- 1 Methods & Experiments**
- 2 Pain measurements**
- 3 Osteogenesis measurements**
- 4 Angiogenesis measurements**

Gait analysis

- Estimate the acute effect of implants in the living body by measuring locomotion



Pain levels can be assessed by measuring locomotion

Gait analysis

- Estimate the acute effect of implants in the living body by measuring locomotion




Pain levels can be assessed by measuring locomotion

- Materials implanted into the right hind (RH) leg



Painful feeling expected to occur in the RH

Gait analysis




• Longer swing time indicates a stronger pain feeling

• ZT1 tends to a lower pain level during recovery which can reduce the sensitivity of mechanical allodynia* induced by the implants

• *Allodynia is extreme sensitivity to touch


- 1 Methods & Experiments**
- 2 Pain measurements**
- 3 Osteogenesis measurements**
- 4 Angiogenesis measurements**

High-resolution microtomography evaluation



- Microarchitecture of the cancellous bone surrounding the **ZT1** is fairly better than CP-Ti and PEEK.

High-resolution microtomography evaluation



• % of bone volume/total volume (BV/TV) significantly higher for ZT1 than PEEK



Higher *in-vivo* osteogenic potential for ZT1 compared with PEEK and CP-Ti

Histological analysis

- Hematoxylineosin (HE) and Masson's trichrome staining

- Hematoxylin stains cell nuclei a purplish blue

- Eosin stains the extracellular matrix and cytoplasm pink

- Masson's trichrome stains collagen, collagen fibers, fibrin, muscles, and erythrocytes.

➡ **Assess the osteointegrations**

Histological analysis


- Hematoxylineosin (HE) and Masson's trichrome staining

- Hematoxylin stains cell nuclei a purplish blue


- Eosin stains the extracellular matrix and cytoplasm pink

- Masson's trichrome stains collagen, collagen fibers, fibrin, muscles, and erythrocytes.

Assess the osteointegrations



Histological analysis



- Quantitative analysis of the bone area ratio (B.Ar/T.Ar)
→ significant differences between the materials


- ZT1 stimulate active bone regeneration and early bone remodeling



Better osteointegration than CP-Ti and PEEK

- 1 Methods & Experiments**
- 2 Pain measurements**
- 3 Osteogenesis measurements**
- 4 Angiogenesis measurements**


Immunofluorescence analysis and imaging of blood vessels



- **Col-I** = osteoblastic marker used to characterize early osteogenesis

- **CD31** = specific endothelial cell surface marker to characterize the angiogenesis

Immunofluorescence analysis and imaging of blood vessels



- **Col-I** = osteoblastic marker used to characterize early osteogenesis


- **CD31** = specific endothelial cell surface marker to characterize the angiogenesis

- **Digital reconstructed radiographs (DRR)** depict the formation of blood vessels around the samples


- More immunofluorescence and more blood vessels are found in the ZT1

Imaging of blood vessels


- **ZT1** shows higher vascular volume than the CP-Ti and the PEEK



It promotes the formation of blood vessels (larger vascular volume postoperatively)




Characterization of the materials



In-vitro measurements



In-vivo measurements



Discussion

ZT1 properties

Biocompatibility


- Wettability
- Chemically stable
- Oxide films for increasing lifetime

- Warning concerning the Cu ions release that can affect blood vessel growth.



What about Biomechanical properties ?


Biomechanical properties



- ZT1 has a more suitable Young's modulus (83 GPa) and a significantly high yield strength (1688 MPa)

Conclusion & Outlooks

Sum-up



Outlooks



Composites Part B: Engineering

Volume 164, 1 May 2019, Pages 800-808



The role of alumina-zirconia loading on the mechanical and biological properties of UHMWPE for biomedical applications

D. Duraccio ^a  , V. Strongone ^{b, c}, G. Malucelli ^b, F. Auriemma ^c, C. De Rosa ^c, F.D. Mussano ^d, T. Genova ^d, M.G. Faga ^a

- Ultra High Molecular Weight Polyethylene (UHMWPE)
- Alumina is used for bone implants (osteoblast interaction)
- Nice frictional and wear behavior
- Nice reliability of tissue ingrowth

Material of choice for hip and knee joint.



**What are the three protagonists
studied ?**

ZT1, PEEK and CP-Ti

**Where was implanted the
material ?**

In the right hind leg

How can we assess the pain level ?

Using gait analysis

How many  were present in
this presentation ?

None
—

Thank you for your attention !

